





Call it a Chair?
SYNCA is much more
than that!



MASSAGE SOFA

PRODUCT SPECIFICATIONS

Product Name	SYNCA Massage Sofa
Color	Beige
Product Weight	70lbs
Japanese Med. Device Cert. Num:	TBD
Weight Limit	250lbs
Wattage	50W
Dimensions	W21.5 x H37 x D39in



JOHNSON FITNESS (M) SDN. BHD. (680946-D)

Lot 557D, Jalan Subang 3, Subang Jaya Industrial Estate, 47610 Subang Jaya, Selangor Darul Ehsan, Malaysia.
TEL: +603-5632 8273 • FAX: +603-5632 7290 • EMAIL: sales@johnsonfitness.com.my
f johnsonfitnessmalaysia • www.johnsonfitness.com.my

THE SPACE OF COMFORT JUST FOR YOU.

SYNCA Massage Sofa is designed from the ground up to give you the most ergonomically comfortable seating angle and position.

The design of the 106-degree angle between the seat and backrest creates a posture so comfortable you won't want to leave.

SYNCA's ambient side lighting and superior massage will set a calming atmosphere every night to help you destress and prepare for a good night's sleep.



SIMPLE TO USE CONTROLLER

Customize your individual massage need to provide all-around solutions or pinpoint massage location for your individual needs.

Location of auto courses

Five auto courses to choose in total: knead and tap the whole body, back focus, waist focus, pelvis and air stretch

Location of manual courses

"Knead & Tap" and "Roller" are two kinds of manual techniques. Combine with 3 body areas: full, upper, and lower body and there are 6 manual massage courses in total.


Location of point courses

The robot can be controlled manually (up and down) to target any sore point precisely.

PROFESSIONAL PHYSIOTHERAPY FOR YOUR PELVIC

Muscles in the pelvic region are an important support structure associated with lower leg and vertebral core muscle movement.

Tense or stiff muscles in pelvic region can cause pain and dysfunction in the lower back and hips.



AT HOME THERAPEUTIC MASSAGE

SYNCA Massage Sofa massages from the buttock to the thigh, carefully tapping, and then deeply kneading to reduce the rigidity and stiffness of the pelvic region. The swaying hip compression from the side air bags further relax the muscles in the pelvic region.



COVERAGE FROM SHOULDER TO THIGH WITH LUMBAR HEAT

The 34-inch SL-Track will massage from the lower neck all the way into the buttocks and hamstring area. Choose full body, back, waist or buttock, or point to focus in specific areas or specific points. It truly is an amazing compact massage chair.

